

CHGD Wenzhou Gaode Medical Instrument Co., Ltd.



Lee: +86 189 6970 2763
Skyla: +86 150 8858 9508
chgd@china-gaode.com



CHGD

Vacuum Blood Collection Tube

Offer Reliable Services For Clinical Examination

Gaode Medical take the concept of

“Safety · Healthy · Environment friendly”

as our Mission.



- **Safety**

Provide human with safe products, to reduce industrial threat to human health.

- **Healthy**

Provide health insurance for families' good times, accurate results keep you away from illness.

- **Environment friendly**

Offer environment friendly products for society to ensure the sustainable development of it.





- Because your work is an important service to people entrusted to your care.
- Because correct analyses results are the basis for far-reaching decisions.
- Because only fully perfected products can guarantee safety and continuous high quality.
- Because more efficiency lessens your burden with routine procedures.





Cap

International Standard Color
Easily Take
Avoid Infection



Stopper

Special silicone
Special material
Avoid block needle



Tube

Special tube silicone
Ensure exact vacuum
Guarantee safe sample



Label

Water-proof treatment
Double label. Bar code
label. OEM label available



Additive


Accurate additive control
Automatic control of amount
A quick-drying technology





- ① CHGD Tube
- ② CHGD Cap
- ③ CHGD Needle- Holder
- ④ Pen-type Blood Collection Needle
- ⑤ Blood Collection Needle with Scalp Vein Set



**How to
Using CHGD Blood Collection system**

1. Ten years' experience of producing **Vacuum Blood Collection Tube**.
2. So far we concentrate on **Vacuum Blood Collection Tube**.
3. Japan **Sekisui** professional engineers' technical guidance.
4. Special frequency **Ultrasonic Wave Cleaning** Technology.
5. Special Inner wall Silicone Technology.
6. Advanced Automatic Production-line.
7. Self-developed **Foil Package** for product protection.





Stopper : Experienced supplier , Various material available to comply with your diagnostic instrument .



Tube : The best supplier in China, help us to control the perfect quality from the start.




Label : Professional manufacturer to satisfy your need. Bar code/double label/OEM label guarantee easy to use.



Additive : Reliable additive supplier help us to provide you pure fine quality, thus to reduce defects during laboratory use.



- Develop purpose:
 - Solve the problem of hemolysis and protein attachment to tubes in clinical examination which caused by unclean inner wall.
- Feature :
 - We use high quality RO-water in ultrasonic wave cleaning equipment to wash the tube.to ensure the clean of tubes, reduce the probability of red cell hemolysis and protein attachment to the tube.



- Develop purpose:
 - Reduce the problems caused by defective raw materials.
- Feature :
 - Inner wall coated with special silicon, This process makes the tubes safe to weak blood cells.
 - Special silicon is stable while surrounded by blood which makes the examination result much more accurate.




- **Develop purpose:**
 - Avoid blood cells hemolysis caused by additives not well spread on the inner wall.
- **Feature :**
 - This production line could control the adding volume of additives, and this machine could spread additives well around the inner wall of tube which makes blood well mixed with additives.



- Develop purpose:
 - Leave much more shelf life to the Sodium Citrate tube.
- Feature :
 - Due to the feature of PET tube, the water will escape after a period. The foil package could stop this event, so the additive could maintain stable for a long shelf time.





➤ Multi-channel you can find us

- Web: www.china-gaode.com/en
- MIC Showroom: <http://gaodemedical.en.made-in-china.com>
- E-mail: chgd@china-gaode.com
- Live chat online: http://www.china-gaode.com/en/service_show.asp?id=162

➤ Technical Support

- Tel: +86 18969702763
- +86 15088589508
- Fax: +86 -577-86851552

➤ Help you contact with other Suppliers in China



CHGD

Your Right Choice

More stable

**We are ready,
what about you?**

More support



More competitive

www.china-gaode.com

Wenzhou Gaode Medical Instrument Co., Ltd.





Product Service

CERTIFICATE

No. Q1N 11 06 66670 003

Holder of Certificate: Wenzhou Gaode Medical Instrument Co., Ltd

Haibin Industrial District
Haibin Street, Longwan
325024 Wenzhou, Zhejiang
PEOPLE'S REPUBLIC OF CHINA

Facility(ies):

Wenzhou Gaode Medical Instrument Co., Ltd
Haibin Industrial District, Haibin Street, Longwan, 325024
Wenzhou, Zhejiang, PEOPLE'S REPUBLIC OF CHINA

Certification Mark:



Scope of Certificate:

Design, Development, Production and
Distribution of Single-use Vacuum
Tubes for Venous Blood Collection

**Applied
Standard(s):**

EN ISO 13485:2003/AC:2009
ISO 13485:2003
Medical Devices – Quality Management Systems –
Requirements for regulatory purposes

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality system which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.:

BJ1196804

Valid from:

2011-06-15

Valid until:

2014-05-31

Date, 2011-06-20

Hans-Heiner Junker



Page 1 of 1

TÜV SÜD Product Service GmbH
Zertifizierstelle
Ridlerstr. 65 - 80339 München
Germany



Akkreditiert durch
Zentralstelle der Länder
für Gesundheitsschutz
bei Arzneimitteln
und Medizinprodukten
ZLG-ZO-999.98.12-46



China

CE Technical Documentation Review Report

(According to IVDD Directive 98/79/EC)

Report No. TFR-BJ-66670 201101

The Manufacturer:

Wenzhou Gaode Medical Instrument Co., Ltd.

Haibin Industrial District, Haibin Street, Longwan, Wenzhou

325024, Zhejiang, P.R.China

Has Established CE Technical Documentation for the Following Products

Single-use Vacuum Tubes for Venous Blood Collection

This report is issued on a voluntary basis. It confirms that the CE technical documentation (no. QR/GDY09-2008, Rev A/0, Date 2010-11-08) of the listed product is complete with regard to the requirements of Annex III section 3 of Directive 98/79/ EC excluding quality system documentation. Detail products covered under the technical documentation are described in this report. This report would be re-assessed before 2012-12-24.

Issue Place:

Medical, Health Service Dept.
TÜV SÜD China, Beijing Branch

Issue Date:

2011-12-24



Any significant changes in design or construction of the product or amendments to the directive or standards applied may render this report invalid. Also, this report is only valid in combination with ISO13485 certificate issued by TÜV SÜD.

Report Number: QIP-ASH113909
Audit Date : 15 Sep., 2011

Expiry Date : 14 Sep., 2012

This report is issued by Focus Technology Co., Ltd. (Made-in-China.com) and the supervising inspectorate (SGS-CSTC Standards Technical Services Co., Ltd.) to confirm that:

Company Name : Wenzhou Gaode Medical Instrument Co., Ltd.
 温州市高德医疗器械有限公司
Showroom : <http://gaodemedical.en.made-in-china.com/>
Address : No. 168, Haining Road, Haibin Industrial Zone, Wenzhou City, Zhejiang, China
Product : Blood Collection Tube, Blood Sampling System, Blood Tube, Medical Consumables, Test Tube, etc.

has been on site audited for the Following Scope of Activity

1. General Information
2. Foreign Trade Capacity
3. Product Research & Development Capacity
4. Quality Management System and Product Certification
5. Working Environment & Energy Saving
6. Photos


General Comments:

Wenzhou Gaode Medical Instrument Co., Ltd. is a manufacturer with 50 employees; it was established in 2003, located in No. 168, Haining Road, Haibin Industrial Zone, Wenzhou City, Zhejiang, China. The workshops occupy an area 3,000 square meters. The company has its own brand "CHGD". They had passed ISO13485:2003 certification, and already had obtained CE certificate for their products.

Sign for and on behalf of
 SGS-CSTC Standards Technical Services Co., Ltd.

This document is issued by the Company under its General Conditions of Service accessible at http://www.sgs.com/terms_and_conditions.htm. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained herein is solely based on visual examination of the subject and is not necessarily representative of the consignment and reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions. If any Client's instructions require a more detailed examination, the responsibility is to its Client and this document does not constitute a guarantee of the accuracy of the information provided. The Company is not liable for any unauthorized alteration, forgery or falsification of the content or appearance of this document in whatever form or by any means. The Company's liability is limited to the amount of the fee paid for the service.
 SGS-CSTC Standards Technical Services Co., Ltd. is a member of the SGS Group. For more information, please contact us at telephone: (86-155)33714431, or email: CN.Docs@sgs.com

